### Produktbeschreibung

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#### Product description

OmniCeram Bulk Fill SDR is a flowable bulk composite with bioceramic fillers for restorations in dentine. OmniCeram Bulk Fill SDR can be light cured in one layer up to 4 mm. OmniCeram Bulk Fill SDR shows optimal flow properties, low polymerization shrinkage, low polymerization stress. The material is ionactive.

#### Indications/Intended use

- Base in cavity class I&II direct restorations
- Liner under direct restorative materials-class II box liner

#### Performance features

The performance features of the product meet the requirements of the intended use.

#### Contraindications

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

#### Patient target group

Persons who are treated during a dental procedure.

#### Intended users

This medical device should only be used by a professionally trained dental practitioner.

#### Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

#### Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

In case of proximal contacts matrix and wedge may be used.

Apply a common adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions. Application with the self-etching bonding OmniCite All'N'One or OmniCite 2.0, which can be used in combination with an etching gel (e.g. Omni Etch), is particularly recommended. For this purpose, the instructions for use for OmniCite All'N'One and OmniCite 2.0 must be observed.

#### Application and Curing of OmniCeram Bulk Fill SDR

Apply OmniCeram Bulk Fill SDR directly into the preparation under slow and steady pressure. Start at the deepest point of the cavity and slowly withdraw the Single Dose Tip with filling the cavity. To prevent air voids. To avoid air bubbles make sure that the Single Dose Tip is always immersed. OmniCeram Bulk Fill SDR shows good flow behavior, a hand instrument may be used if necessary.

When used als bulk fill base material, most cavities can be filled in one increment (up to 4 mm). Leave at least 2 mm for the occlusal filling composite.

In case of deeper cavities apply OmniCeram Bulk Fill SDR in increments of 4 mm and cure each layer separately.

It is also possible to apply OmniCeram Bulk Fill SDR as liner in a thin layer onto exposed dentine areas.

By using a polymerization unit (wavelength range 400-500 nm) with a **light intensity of at least 1000mW/cm<sup>2</sup>**, cure each increment for **20 seconds**.

Hold the light emission window as close as possible to the filling material.

Complete the restoration with a suitable methacrylate based light cure composite according to the corresponding instructions.

#### Finishing/Polishing

If finishing and/or polishing of the cured OmniCeram Bulk Fill SDR is necessary, use traditional rotary instruments according to manufacturer instructions.

#### Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

#### Additional Notes/Warnings

- If the material is applied directly from a Single Dose Tip into the mouth, the Single Dose Tip is for single use only due to hygienic reasons.
- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

#### Composition

Dimethacrylates, dental glass, silicon dioxide, photo initiators.

#### Disposal

Disposal of the product according to local authority regulations.

#### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

#### Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

#### Warranty

SanaPro Dental GmbH warrants this product will be free from defects in material and manufacture. SanaPro Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and SanaPro Dental GmbH's sole obligation shall be repair or replacement of the SanaPro Dental GmbH product.

#### Limitation of Liability

Except where prohibited by law, SanaPro Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

SanaPro Dental GmbH  
Stresemannstrasse 46  
D-27570 Bremerhaven  
Germany

OMNIDENTAL-  
Handels-gesellschaft mbH  
Gutenbergring 5  
D-63110 Rodgau