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**OmniCeram multi shade** is a light cure radiopaque composite for anterior and posterior restorations. OmniCeram multi shade can be used for all natural tooth shades due to its continuous color matching properties. Beside its permanent gloss and toothlike opalescence OmniCeram multi shade shows excellent physical properties and good handling features. High gloss polishability demands only little effort. OmniCeram multi shade contains dimethacrylate based resins and inorganic fillers < 0.2 µm. The total filler content is 80 % by weight and 75 % by volume. Radiopacity of OmniCeram multi shade is 200 % aluminium (radiopacity of enamel: ~200 % aluminium, of dentine: ~ 100 % aluminium).

**OmniCeram multi shade dentin** is a light cure composite that is used in combination with OmniCeram multi shade to mask discolored teeth. OmniCeram multi shade dentin contains dimethacrylate based resins and inorganic fillers < 0.2 µm. The total filler content is 80 % by weight and 76 % by volume. Radiopacity of OmniCeram multi shade dentin is 200 % aluminium (radiopacity of enamel: ~200 % aluminium, of dentine: ~ 100 % aluminium).

The material properties of OmniCeram multi shade and OmniCeram multi shade dentin meet the requirements of **ISO 4049**, type 1, class 2, group 1.

#### Indications OmniCeram multi shade

- Restorations of all cavity classes
- Direct veneers
- Diastema closure
- Repair of ceramic restorations (e.g. veneers, crowns)

#### Indications OmniCeram multi shade dentin

- Masking of the oral cavity in case of extended class III and IV cavities
- Masking of tooth discoloration

#### Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with pulp capping material.

#### Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

#### Preparation of Tooth

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

#### Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. For pulp protection areas close to the pulp should be covered with a pulp capping material. Refer to the corresponding user instructions. Apply a common adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions. Recommended adhesives are OmniCite All'N'One, a self etching adhesive or OmniCite 2.0, an adhesive for the Total-Etch-Technique in combination with OmniEtch.

#### Application and Curing of OmniCeram multi shade

In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure flowable composite (e.g. OmniCeram evolution flow 2.0), followed by the application of OmniCeram multi shade.

In case of small cavities OmniCeram multi shade can be applied directly.

#### Note:

In case of extended class III and class IV cavities use OmniCeram multi shade dentin palatinal to ensure that the restoration may not appear dark. Use OmniCeram multi shade dentin also as a first layer for class I, II and V cavities of dark or discolored teeth.

In case of repair of ceramic or composite restorations use a light cure opaquer as first layer in case of metal or zirconia surfaces.

Apply OmniCeram multi shade dentin in **layers of max. 1 mm**. Light cure each increment separately.

**By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm<sup>2</sup>, cure each increment separately for 20 seconds.**

Hold the light emission window as close as possible to the filling material.

Apply OmniCeram multi shade in **layers of max. 2 mm**. Light cure each increment separately. OmniCeram multi shade is easy to model. Place it by using a suitable instrument. Matrix strips may be used.

**By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm<sup>2</sup>, cure each increment separately for 20 seconds.**

Hold the light emission window as close as possible to the filling material.

#### Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs.

Polish to high gloss with polishing discs or polishing brushes.

#### Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

#### Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

#### Warranty

SanaPro Dental GmbH warrants this product will be free from defects in material and manufacture. SanaPro Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and SanaPro Dental GmbH's sole obligation shall be repair or replacement of the SanaPro Dental GmbH product.

#### Limitation of Liability

Except where prohibited by law, SanaPro Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

**Keep away from children!  
For dental use only!**

 **Manufacturer:**  
SanaPro Dental GmbH  
D-27570 Bremerhaven  
Germany

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