## Omnidrill

Universal Bohrebad, aldehydfrei

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### GB Instructions for use

#### Type of application:

For the cleaning and prophylactic disinfection of drills, cutters, diamond ringrinders and other rotating instruments. Not suitable for the terminal disinfection of semi-critical and critical medical devices!

#### Notes on application:

Always use Omnidrill undiluted! After use, put the instruments straight into the ready-for-use solution, all the parts to be disinfected have to be fully covered. After disinfection the instruments should be rinsed, ideally with demineralised water, and dried thoroughly. Do not mix with other disinfectants or cleaners, not suitable for aluminium. Depending on the degree of contamination, the working solution may need to be replaced (RKI recommendation: at least once per working day). The working solution can be disposed of in normal drainage. It is recommended to wear safety gloves and eye protection.

#### Scope of effectiveness:

Bactericidal (incl. MRSA) as per EN 13727, EN 14561, VAH; levurocidal (Candida albicans) as per EN 13624, EN 14562, VAH; tuberculocidal (M. terrae) as per EN 14348; limited virucidal effect (enveloped viruses such as HBV, HIV, HCV) as per DVV/RKI..

#### Application times:

Bactericidal (incl. MRSA) according to EN 13727, EN 14561, VAH, dirty conditions): Conc.: 5 min.  
Levurocidal (Candida albicans (EN 13624, EN 14562, VAH, dirty conditions): Conc.: 5 min.  
Tuberculocidal (M. terrae) as per EN 14348, EN 14563, clean conditions: Conc.: 5 min.  
Limited virucidal effect (Enveloped viruses (e.g. HBV, HCV and HIV) as per RKI/DVV (dirty conditions): Conc.: 5 min.

#### 100g contain:

100g contains: 0.396g didecyldimethylammonium chloride, 0.2475g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamin.

#### Safety guidelines:

H412 Harmful to aquatic life with long lasting effects. Avoid release to the environment. Dispose of the contents/container in accordance with local/regional/ national/international regulations.

Directive 648/2004/EC: Contains < 5% non-ionic surfactants, <5% cationic surfactants, <5% EDTA salts, fragrance.

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PRISMAN GmbH